PRODUCT BROCHURE

## Digital Angel Radar



CARING FURTHER
ALPHALIFECARE.COM.AU



## Falls Detection Technology Comparison



	Alpha Lifecare Radar Solution	Other Radar	Legacy Device
Falls Detection	~	<b>✓</b>	<b>✓</b>
Falls Prevention	~	Limited	No
Falls Analytics	~	Limited	No
ADL Analytics - daily routines	~	No	No
Gait Analytics - falls prediction, detection of changes	✓	No	No
Area Coverage	40 Sqm	15 Sqm	5 Sqm
See Through the Walls	<b>✓</b>	No	No
ML and Al Analytics - near Zero false alerts rate	✓	No	No
Smart Notifications	<b>✓</b>	<b>✓</b>	No
Respiratory Rate	<b>✓</b>	No	No
Heart Rate Monitoring	~	No	No
Sleep Review Analytics	~	No	No
Secondary Alerts - long laying on the floor, seating in bed	✓	No	No
Full Dashboard with Analytics	<b>✓</b>	Third Party	Third Party
iOS and Android App	~	Third Party	<b>~</b>
Integrate with Nurse Call Solution	<b>✓</b>	Limited	No
Visitor Numbers	<b>✓</b>	Limited	No
Chair Alerts	<b>✓</b>	No	No

### Less Sensors, More Care



### Alpha Lifecare Radar Solution



- One device covers the entire room/bathroom
- · Discreet recessed model
- Ability to monitor activity levels
- High accuracy levels
- Ability to integrate with other solutions

### Other Systems



- Increased quantity of sensors per room
- Sensors placed around the room
- Limited ability to monitor activity levels
- High rate of false positives

### Overview



#### Alpha Lifecare Radar Solution

- Movement
- Gait
- Precise location
- Respiratory
- Heart rate
- Sleep analysis
- Falls detection
- Entire room

### Alpha Lifecare Radar Solution Devices

- Can detect falls from Floorline bed
- Detection is not affected by air mattress operation (incl. respiratory & heart rate)
- No interference from ceiling fans, ceiling lift, or A/C units
- 100% Non Invasive (No Cameras or wearables)

### Other Radar

- Movement
- Falls detection
- Separate devices must be installed in each room and bathroom
- Wall Placement
- Battery operated in some cases


### Alpha Lifecare Radar Solution Devices

- Can detect falls from Floorline bed
- Detection is not affected by air mattress operation (incl. respiratory & heart rate)
- No interference from ceiling fans, ceiling lift, or A/C units
- 100% Non Invasive (No Cameras or wearables)

### Respiration Monitoring



## Accurate Respiration Parameters Extraction Advanced and patented analysis approach



#### **Parameters**

- Respiration Rate
- Inhale/Exhale Time
- Stop Time


## Abnormal Cases Detection

- Hyperventilation
- Apnea

### Smart Dashboard



Included in your subscription and can be integrated with EHR and nurse call via API



### Validated Proprietary Technology



With strong IP portfolio that includes 6 issued patents, Alpha Lifecare Radar Solution leads the category with unique validated proprietary technology.

Alpha Lifecare developed industry leading RF based system that enables monitoring 5 body Indicators constantly measured followed by AI & ML Engines.







**Posture** 




Motion



Respiration



**Heart Rate** 



## Alpha Cloud Includes:

- Access to healthcare professional
- API for other operators' NCS
- Remote SW update and new services
- Big Data analytics

# How It Works



#### 4 Body Indicators Constantly Measured:







Posture



Motion



Respiration

#### The Detection Process:

RF Interferometer (RADAR)

Human States Classifier

Cognitive Situation Analysis (AI, ML)

:

WiFi Communication

Cloud / Dashboard

•

**Alert Given** 

## We Look Forward To Working With You



Alpha Lifecare  $\Lambda$ 

1300 930 930 info@alphalifecare.com.au